

ProteoSpin™ Detergent Clean-Up Micro Kit

The ProteoSpin™ Detergent Clean-up Micro Kit provides a fast and simple procedure for the effective removal of SDS, Triton® X-100 and other detergents from small protein samples. The kit employs spin-column chromatography using Norgen's patented protein resin as an ion exchanger. Each column is capable of processing up to 50 µg of either acidic or basic protein. The ProteoSpin™ Detergent Clean-up Micro Kit provides a convenient method for removing detergents from small protein samples before running downstream applications including mass spectrometry, SDS-PAGE, proteolytic digestion, isoelectric focusing, x-ray crystallography and NMR spectroscopy.



Kit Specifications			
Maximum Protein Input	50 µg	% Detergent Removal	> 95%
Protein Recovery	80 - 95%	Minimum Elution Volume	30 µL
Column Binding Capacity	50 µg	Time to Process 12 Samples	20 minutes

ProteoSpin™ Detergent Clean-up Micro Kit Benefits

Fast and easy processing	This easy-to-use procedure allows processing of up to 12 samples in only 20 minutes
No detergent carryover	Each sample is processed individually using its own spin column, resulting in no carryover or column bleed from sample to sample.
Useable with a variety of detergents	The kit can be used to remove detergents including SDS, Triton® X-100, CHAPS, NP-40, and Tween 20.
Effective detergent removal	Removes more than 95% of detergents for many protein samples, allowing tryptic digestion.
High recovery	Protein recoveries of 80-95% (using initial detergent concentrations of 0.1%) for many proteins (Figure 1).
Suitable for downstream applications	Final elution suitable for: <ul style="list-style-type: none"> • Mass spectrometry • SDS-PAGE • Proteolytic digestion • Isoelectric focusing • X-ray crystallography • NMR spectroscopy

ProteoSpin™ Detergent Clean-up Micro Kit

Sample	A ₆₀₀	% SDS Remaining	% SDS Removal	mg SDS Removed
Control (water)	0	0	0	0
0.1% SDS	0.142			
ProteoSpin sample 1	0.024	0.017%	99.98%	2.457
ProteoSpin sample 2	0.030	0.02%	99.98%	2.458

Table 1. Efficient Removal of SDS

Removal of SDS with the ProteoSpin™ Detergent Clean-up Micro Kit is greater than 99%.

Detergent Clean-up Micro Kit Contents

1. Column Activation and Wash Buffer for Acidic Proteins
2. pH Binding Buffer for Acidic Proteins
3. Column Activation and Wash Buffer for Basic Proteins
4. pH Binding Buffer for Basic Proteins
5. Elution Buffer
6. Neutralizer
7. Micro Spin Columns
8. Elution Tubes
9. Application Manual

Shipping Conditions

The ProteoSpin™ Detergent Clean-up Micro Kit is shipped at room temperature.




Figure 1. High Protein Recovery in the Presence of SDS

Solutions containing 50 µg BSA and 0-1.0% SDS were processed using the ProteoSpin™ Detergent Clean-up Micro Kit to remove SDS. High protein recovery was obtained even with original SDS concentration as high as 1%.

Customer-Supplied Reagents and Equipment

- Benchtop microcentrifuge
- Micropipettors
- pH indicator paper
- Other optional elution buffers
- Isopropanol
- Sterile, deionized water or Milli-Q® water

Storage Conditions

Unopened solutions should be stored at room temperature. Once opened, the solutions should be stored at 4°C when not in use, except for the Basic and Acidic Binding Buffers, which should be stored at room temperature.

Cat #	Description	Quantity
10200	ProteoSpin™ Detergent Clean-up Micro Kit	25 preps
10500	ProteoSpin™ Detergent Clean-Up Micro Kit	50 preps