

Saliva DNA Collection, Preservation and Isolation Kit – 50 Individual Devices
Product Insert

Product #: RU35700

Norgen’s Saliva DNA Collection, Preservation and Isolation Kit provides a safe and rapid all-in-one procedure for the collection, preservation and isolation of saliva DNA at ambient temperatures.

INTENDED USE

Norgen’s Saliva DNA Collection, Preservation and Isolation Kit is an all-in-one solution designed for 1) simple and non-invasive saliva collection; 2) preservation of DNA in saliva samples at ambient temperature; and 3) isolation of high quality DNA within a laboratory setting. The Saliva DNA Collection, Preservation and Isolation Kit contains 50 Individual Saliva DNA Collection and Preservation Devices, as well as the required reagents for the subsequent laboratory isolation of the saliva DNA from the preserved samples. Each of the 50 individual Saliva DNA Collection and Preservation Devices consists of 3 components: (1) Saliva Collection Funnel and Collection Tube, (2) Collection Tube Cap, and (3) Norgen’s Saliva DNA Preservative contained within a sealed squeezable ampoule. Saliva samples are collected by spitting inside the Collection Funnel which has been assembled with the Collection Tube. After collecting the required volume of saliva the Collection Funnel is removed and the contents of the Preservative Ampoule are then added and mixed with the collected saliva. The Saliva Collection Tube is subsequently sent to the laboratory for DNA isolation and analysis. Each of Norgen’s Collection Tubes is labeled with a unique serial number that can be used for secure and anonymous tracking of the sample. The saliva DNA in preserved samples is stable for more than 2 years at room temperature. This kit is ideal for collecting, preserving and isolating DNA samples for epidemiological and population studies.

SALIVA DNA PRESERVATIVE

Norgen’s Saliva DNA Preservative is an aqueous storage buffer designed for rapid cellular lysis and subsequent preservation of DNA from fresh specimens. The buffer prevents the growth of Gram-negative and Gram-positive bacteria and fungi, and also inactivates viruses allowing the resulting non-infectious samples to be handled and shipped safely. In addition, the buffer eliminates the need to immediately process or freeze samples and allows the samples to be shipped to centralized testing facilities at ambient temperature. The components of the buffer allow samples to be stored for more than 2 years without any detectable DNA degradation.

NUCLEIC ACID ISOLATION FROM PRESERVATIVE

Saliva DNA is isolated from the preserved saliva samples via alcohol precipitation using the provided reagents. Saliva DNA can also be isolated from the preserved saliva samples using Norgen’s spin-column based Saliva DNA Isolation Kit (Cat# RU45400).

Saliva DNA Collection, Preservation and Isolation Kit Contents:

Component	Contents
Individual Saliva DNA Collection and Preservation Device	50
Proteinase K	1 vial
Product Insert	1

Individual Saliva DNA Collection and Preservation Device Contents:

Component	Contents
Saliva Collection Funnel and Collection Tube	1
Collection Tube Cap	1
Preservative Ampoule	1
Donor Instructions	1

Kit Specifications	
Volume of Saliva Collected	2 mL
Volume of Saliva-Preservative Mix	4 mL
Preservation Temperature	Room temperature
Preservation Time	Over 2 years at room temperature
Time to Complete DNA Isolation	45 minutes plus incubation time
Average Yield from 0.5 mL*	20 µg
Average OD 260/280	1.7

* Average DNA yield will vary depending on the health status of the donor

MATERIALS REQUIRED BUT NOT SUPPLIED

For DNA Isolation

- Benchtop microcentrifuge
- Micropipettors
- 2 mL microcentrifuge tubes
- 55°C incubator
- 70% ethanol
- Isopropanol
- TE Buffer (10 mM Tris pH 8.0, 1 mM EDTA)

SHELF LIFE AND HANDLING

- The **Collection Device** should be kept tightly sealed and stored at room temperature for up to 3 years without any reduction in kit performance (The collection due date is written on the device label).
- Once collected, saliva is stable for more than 2 years when kept tightly sealed and stored at room temperature.
- The kit contains ready-to-use **Proteinase K** which is dissolved in a specially prepared storage buffer. The Proteinase K should be stored at -20°C for up to 2 years after the “Collect Before Date” (with minimum freeze and thaw).
- The Collection Tube, the Collection Funnel and the Device Container of each Individual Saliva DNA Collection and Preservation Device are recyclable.

QUALITY CONTROL

In accordance with Norgen’s Quality Management System, each lot of Norgen’s Saliva DNA Collection, Preservation and Isolation Kit is tested against predetermined specifications to ensure consistent product quality.

Disclaimers and Satisfaction Guarantee

This kit is designed for research purposes only. It is not intended for diagnostic use.

NORGEN BIOTEK CORPORATION guarantees the performance of all products in the manner described in our product manual. The customer must determine the suitability of the product for its particular use. We reserve the right to change, alter, or modify any product to enhance its performance and design.

WARNINGS AND PRECAUTIONS – Saliva Donor

1. If any of the Saliva DNA Preservative comes into contact with skin or eyes, wash thoroughly with water. Do not ingest the Saliva DNA Preservative.
2. The collection device contains a choking hazard (Collection Tube Cap). Keep away from children.

WARNINGS AND PRECAUTIONS – DNA Isolation from Preserved Saliva Samples

- 1 Follow universal precautions. All specimens should be considered as potentially infectious and handled accordingly.
- 2 Wear personal protective equipment, including gloves and lab coats when handling kit reagents. Wash hands thoroughly when finished performing the procedures.
- 3 Do not smoke, drink or eat in areas where kit reagents and/or human specimens are being used.
- 4 Dispose of unused kit reagents and human specimens according to local, provincial or federal regulations.
- 5 As contamination of specimens or reagents can produce erroneous results, it is essential to use aseptic techniques.
- 6 Only use the protocol provided in this insert. Alterations to the protocol and deviations from the times and temperatures specified may lead to erroneous results.
- 7 For more information, please consult the appropriate Material Safety Data Sheets (MSDSs). These are available as convenient PDF files online at www.norgenbiotek.com.

Procedures

A. SAMPLE COLLECTION AND PRESERVATION PROCEDURES

- 1 Do not eat, chew gum, drink or smoke for 30 minutes prior to saliva collection.
- 2 Rinse your mouth with water and wait for 10 minutes.
- 3 Spit into the **Collection Funnel** until the amount of liquid saliva (not including bubbles) in the **Collection Tube** reaches the 2 mL line. Discard / recycle the Collection Funnel.
- 4 Twist and remove the tip of the **Preservative Ampoule**.
- 5 Squeeze the **Preservative Ampoule** to empty the contents into the **Collection Tube**. Screw on the cap of the **Collection Tube** to close it tightly.
- 6 Shake the **Collection Tube** well for 10 seconds to mix saliva and preservative.
- 7 The saliva sample is now preserved and ready for storage, shipping or processing.
Note: Each of Norgen's Collection Tubes is labeled with a unique serial number that can be used for secure and anonymous tracking of the sample.
- 8 Dispose of the **Preservation Ampoule** and tip in the garbage. The plastic device container can be recycled.

B. INSTRUCTIONS FOR STORAGE OF PRESERVED SALIVA SAMPLES

Preserved saliva samples can be stored at room temperature for more than 2 years without significant loss of DNA quality.

C. ISOLATION OF SALIVA DNA FROM 0.5 mL OF PRESERVED SALIVA SAMPLES

Note: This procedure outlines the isolation of DNA from 0.5 mL of the preserved saliva samples. The procedure for the isolation of DNA from 4 mL of preserved saliva samples is in Section D. Please note that procedures for the isolation of DNA using high throughput 96-well formats are available on Norgen's website at www.norgenbiotek.com


1. Incubate the preserved saliva sample in the **Collection Tube** at 55°C for 30 minutes in a water incubator or 1 hour in an air incubator.
2. Mix the preserved saliva sample by inversion and gentle shaking for a few seconds.
3. Transfer 0.5 mL of preserved saliva sample to a 2 mL centrifuge tube (not provided).
4. Add 10 µL of Proteinase K, mix by vortexing for a few seconds and incubate the sample at 55°C for 1 hour.
5. Add an equal volume of room temperature isopropanol to the sample. Mix gently by inversion 10 times.
6. Centrifuge at room temperature for 5 minutes at maximum speed (14,000 xg).
7. Carefully remove and discard the supernatant, taking care not to disturb the DNA pellet. Gently place the tube briefly upside down on a paper towel to remove residual isopropanol.
8. Carefully add 500 µL of 70% ethanol. Gently swirl and let stand at room temperature for 1 minute.
9. Centrifuge at room temperature for 1 minute at maximum speed (14,000 x g), then carefully remove and discard the ethanol without disturbing the pellet.
10. Gently place the open tube upside down on paper towel for 10-15 minutes to remove the excess amount of 70% EtOH and to air dry the DNA pellet.
11. Add 100 µL of TE Buffer (not provided). Vortex for 30 seconds and incubate at 55°C for 5 minutes to rehydrate the DNA pellet. Ensure complete rehydration of the DNA (pellet and smear on the side of the tube) prior to any subsequent step.

Note: Large amounts of high molecular DNA can be very slow to rehydrate (dissolve) completely. Incomplete rehydration of the DNA will affect estimating the DNA concentration and can result in failure of downstream applications such as SNP, digestion and PCR.

12. Centrifuge the tube at 14,000 x g for 1 minute to pellet any insoluble material.
13. Transfer the clear liquid into a clean tube taking care not to disturb the pellet.
14. The purified DNA sample may be stored at 4°C up to 2 months. It is recommended that samples be placed at -20°C for long term storage.

Flowchart

Procedure for Isolating Saliva DNA from 0.5 mL of Preserved Saliva Samples



D. ISOLATION OF SALIVA DNA FROM 4 mL OF PRESERVED SALIVA SAMPLES


Note: This procedure outlines the isolation of DNA from the entire volume preserved saliva samples (4 mL). Please note that procedures for the isolation of DNA using high throughput 96-well formats are available on Norgen's website at www.norgenbiotek.com

1. Incubate the preserved saliva sample in the **Collection Tube** sample at 55 °C for 30 minutes in a water incubator or 1 hour in an air incubator.
2. Mix the preserved saliva sample by inversion and gentle shaking for a few seconds.
3. Transfer 4 mL of preserved saliva sample to a 15 mL centrifuge tube (not provided).
4. Add 80 µL of Proteinase K, mix by vortexing for a few seconds and incubate the sample at 55 °C for 1 hour.
5. Add an equal volume of room temperature isopropanol to the sample. Mix gently by inversion 10 times.
6. Centrifuge at room temperature for 10 minutes at a minimum of 3,500 xg.
7. Carefully remove and discard the supernatant, taking care not to disturb the DNA pellet. Gently place the tube briefly upside down on a paper towel to remove residual isopropanol.
8. Carefully add 1 mL of 70% ethanol. Let stand at room temperature for 1 minute.
9. Centrifuge at room temperature for 5 minutes at a minimum of 3,500 x g, then carefully remove and discard the 70% ethanol without disturbing the pellet.
10. Gently place the open tube upside down on paper towel for 10-15 minutes to remove the excess amount of 70% EtOH and to air dry the DNA pellet.
11. Rehydrate the DNA pellet in 500 µL of TE Buffer (not provided). Pipette or vortex to dissolve the pellet. Ensure complete rehydration of the DNA (pellet and smear on the side of the tube) prior to any subsequent step. Incubate at 55 °C for 30 minutes with occasional vortexing

Note: Large amounts of high molecular DNA can be very slow to rehydrate (dissolve) completely. Incomplete rehydration of the DNA will affect estimating the DNA concentration and can result in failure of downstream applications such as SNP, digestion and PCR.
12. Transfer all the liquid to a 1.5 mL microcentrifuge tube and centrifuge at 14,000 x g for 1 minute to pellet any insoluble material.
13. Transfer the clear liquid into a clean tube taking care not to disturb the pellet.
14. The purified DNA sample may be stored at 4 °C up to 2 months. It is recommended that samples be placed at -20 °C for long term storage.

Flowchart

Procedure for Isolating Saliva DNA from 4 mL of Preserved Saliva Samples



Related Products	Product #
Shipping Accessories – 50 Mailers	38300
Saliva DNA Isolation Kit	RU45400
Saliva DNA Collection and Preservation Devices (50)	RU49000

Technical Assistance

NORGEN's Technical Service Department is staffed by experienced scientists with extensive practical and theoretical expertise in sample and assay technologies and the use of NORGEN products. If you have any questions or experience any difficulties regarding Norgen's Saliva DNA Collection, Preservation and Isolation Kits or NORGEN products in general, please do not hesitate to contact us.

NORGEN customers are a valuable source of information regarding advanced or specialized uses of our products. This information is helpful to other scientists as well as to the researchers at NORGEN. We therefore encourage you to contact us if you have any suggestions about product performance or new applications and techniques.

For technical assistance and more information, please contact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362 or call one of the NORGEN local distributors (www.norgenbiotek.com) or through email at techsupport@norgenbiotek.com.

3430 Schmon Parkway, Thorold, ON Canada L2V 4Y6
 Phone: (905) 227-8848
 Fax: (905) 227-1061
 Toll Free in North America: 1-866-667-4362