



NEW MES and MOPS Running Buffers

MES-SDS-Running Buffer (20X)

Composition

MES, 1M; Tris Base, 1M; SDS, 2%; EDTA, 20mM, pH 7.3 ± 0.1

Method

Dilute 20X with de-ionised water before use; slight colouration with time will not affect the quality

MOPS-SDS Running Buffer (20X)

Composition

MOPS, 1M; Tris Base, 1M; SDS, 2%; EDTA, 20mM, pH 7.7 ± 0.1

Method

Dilute 20X before use; slight colouration with time will not affect the quality

Ordering Information

| Catalogue Number | Description | List Price |
|------------------|--|------------|
| B9-0092 | MOPS SDS Running Buffer (20X), 500ml | 41.55 |
| B9-0094 | MOPS SDS Running Buffer (20X), 1 Litre | 79.00 |
| B9-0096 | MES SDS Running Buffer (20X), 500ml | 41.55 |
| B9-0098 | MES SDS Running Buffer (20X), 1 Litre | 79.00 |