



## NEW MES and MOPS Running Buffers

### MES-SDS-Running Buffer (20X)

#### Composition

MES, 1M; Tris Base, 1M; SDS, 2%; EDTA, 20mM, pH 7.3 ± 0.1

#### Method

Dilute 20X with de-ionised water before use; slight colouration with time will not affect the quality

### MOPS-SDS Running Buffer (20X)

#### Composition

MOPS, 1M; Tris Base, 1M; SDS, 2%; EDTA, 20mM, pH 7.7 ± 0.1

#### Method

Dilute 20X before use; slight colouration with time will not affect the quality

## Ordering Information

Catalogue Number	Description
B9-0092	MOPS SDS Running Buffer (20X), 500ml
B9-0094	MOPS SDS Running Buffer (20X), 1 Litre
B9-0096	MES SDS Running Buffer (20X), 500ml
B9-0098	MES SDS Running Buffer (20X), 1 Litre