



NEW MES and MOPS Running Buffers

MES-SDS-Running Buffer (20X)

Composition

MES, 1M; Tris Base, 1M; SDS, 2%; EDTA, 20mM, pH 7.3 ± 0.1

Method

Dilute 20X with de-ionised water before use; slight colouration with time will not affect the quality

MOPS-SDS Running Buffer (20X)

Composition

MOPS, 1M; Tris Base, 1M; SDS, 2%; EDTA, 20mM, pH 7.7 ± 0.1

Method

Dilute 20X before use; slight colouration with time will not affect the quality

Ordering Information

Catalogue Number	Description	List Price
B9-0092	MOPS SDS Running Buffer (20X), 500ml	41.55
B9-0094	MOPS SDS Running Buffer (20X), 1 Litre	79.00
B9-0096	MES SDS Running Buffer (20X), 500ml	41.55
B9-0098	MES SDS Running Buffer (20X), 1 Litre	79.00