



***TransIT-HeLaMONSTER*® Transfection Kit**

A high efficiency, low toxicity, DNA transfection reagent optimized for HeLa cells

- **Unique Formulation**—Each kit contains a base and booster reagent to maximize transfection performance in HeLa cells.
- **Low Cellular Toxicity**—Maintain cell density and reduce experimental biases.
- **High Efficiency Delivery**—Achieve 50-60% transfection efficiency to ensure experimental success.
- **Deliver Single or Multiple Plasmids**—Ideal for specialized applications, such as gene expression and promoter analysis.

Data

Figure 1. The *TransIT-HeLaMONSTER* Kit Achieves High Transfection Efficiency.

HeLa cells transfected with pEGFP using *TransIT-HeLaMONSTER* Kit in complete media for 24 hours.

